



P.O. Box 7307, North Brunswick, NJ 08902-7307  
Tel: 732-819-0381 Fax: 732-777-5129 E: [sanni@globepharma.com](mailto:sanni@globepharma.com)  
[www.globepharma.com](http://www.globepharma.com), [www.maxiblend.net](http://www.maxiblend.net)

## CAPSULE POLISHER MODEL GPCP-350-WD

GlobePharma's Model GPCP-350-WD is a versatile polisher, which can be used for both capsule polishing and tablet de-dusting. It has a rotating brush with built-in GMP features that facilitate use, cleaning and maintenance. It can be washed down with water.



### FEATURES:

- All product contact parts are made of stainless steel
- Rate of output of 200,000 plus capsules per hour
- Variable speed operation
- Snap and lift assembly and disassembly
- Sealed ball bearings
- Seamless Nylon polishing net which slides easily inside the perforated tube
- Compressed air is provided along the entire length of the finishing chamber
- Coil Spring, positioned inside the chamber ensures continuous, forward movement of heavy and large product
- All replacement parts including brush, net, and shaft can be replaced individually without need to replace the whole assembly
- Easy cleaning and maintenance
- All the electrical, pneumatic and the replacement parts are US made.
- Can be washed down with water.

## **SPECIFICATIONS:**

### **Electrical:**

NEMA 4 DC drive is mounted to a ¼ HP NEMA 4 90 VDC, 1800 rpm motor.

AC Volts:	120, 50/60 Hz
AC Amps:	4.0
DC Volts:	0-90
DC Amps:	2.5

### **Physical Data:**

Width:	24"
Height:	40"
Depth:	30"
Discharge height (adjustable)	27"-30"
Inlet height (adjustable)	33"-38"
Weight net:	55 lb
Weight shipped:	120 lb

### **Materials of Construction:**

All the product contact parts are made of stainless steel and the rest of the machine is also substantially made of stainless steel. The base of the machine is made of polished chrome. The screen is made of FDA approvable Nylon.